

Integrin beta 3(N-terminus) Monoclonal Antibody

Catalog No: #27094

Package Size: #27094-1 50ul #27094-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Integrin beta 3(N-terminus) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3C12-G2-E4-H9
Isotype	IgG1
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu D Rb Rt Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human ITGB3 protein fragments expressed in E.coli.
Target Name	Integrin beta3
Other Names	GT, CD61, GP3A, BDPLT2, GPIIIa
Accession No.	Swiss-Prot:P53041 Gene ID:5536
Uniprot	P53041
GeneID	5536;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol

Application Details

Predicted MW: 120kd

Western blotting: 1:1000

Images



Western blot detection of ITGB3 in Leukocyte whole cell lysate(1:1000 diluted).
 Predicted band size:120KDa.
 Observed band size:120KDa.

Background

Integrin alpha-V/beta-3 is a receptor for cytotactin, fibronectin, laminin, matrix metalloproteinase-2, osteopontin, osteomodulin, prothrombin,

thrombospondin, vitronectin and von Willebrand factor. Integrin alpha-IIb/beta-3 is a receptor for fibronectin.

Note: This product is for in vitro research use only