

## NUP98 Monoclonal Antibody

Catalog No: #27111

Package Size: #27111-1 50ul #27111-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	NUP98 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3B8-D7-H10
Isotype	IgG1
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu D Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human NUP98 protein fragments expressed in E.coli.
Target Name	NUP98
Other Names	ADIR2, NUP96, NUP196
Accession No.	Swiss-Prot:P31939Gene ID:471
Uniprot	P31939
GeneID	471;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

## Application Details

Predicted MW: 98kd


Western blotting: 1:1000

## Images



Western blot detection of NUP98 in HeLa

Immunoprecipitation analysis of HeLa cell lysates using NUP98 antibody



## Background

Nup98 and Nup96 play a role in the bidirectional transport across the nucleoporin complex (NPC). The repeat domain in Nup98 has a direct role in the transport.

Note: This product is for in vitro research use only