

LRRFIP1 Monoclonal Antibody

Catalog No: #27113

Package Size: #27113-1 50ul #27113-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	LRRFIP1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4F3-D8-C9
Isotype	IgG1
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human LRRFIP1 protein fragments expressed in E.coli.
Target Name	LRRFIP1
Other Names	GCF2, TRIP, FLAP1, GCF-2, FLAP-1, HUF1-1, FLIAP1
Accession No.	Swiss-Prot:O43852Gene ID:813
Uniprot	O43852
GeneID	813;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol

Application Details

Predicted MW: 160kd

Western blotting: 1:1000

Images



Western blot detection of LRRFIP1 in Raji, C2C12, MCF7 & Jurkat cell lysates using LRRFIP1 antibody (1:1000 diluted).
 Predicted band size: 89KDa.
 Observed band size: 160KDa

Background

Transcriptional repressor which preferentially binds to the GC-rich consensus sequence (5'-AGCCCCGGCG-3') and may regulate expression of

TNF, EGFR and PDGFA. May control smooth muscle cells proliferation following artery injury through PDGFA repression.

Note: This product is for in vitro research use only