

ATIC Monoclonal Antibody

Catalog No: #27115

Package Size: #27115-1 50ul #27115-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	ATIC Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3H12-C9-H9
Isotype	IgG1
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu Rb
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human ATIC protein fragments expressed in E.coli.
Target Name	ATIC
Other Names	PURH, AICAR, AICARFT, IMPCHASE
Accession No.	Swiss-Prot:P06239Gene ID:3932
Uniprot	P06239
GeneID	3932;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol

Application Details

Predicted MW: 64kd

Western blotting: 1:1000

Images



Western blot detection of ATIC in various tissue and cell lysates using ATIC antibody (1:1000 diluted).
 Predicted band size:64KDa.
 Observed band size:64KDa

Background

Bifunctional enzyme that catalyzes 2 steps in purine biosynthesis.

Note: This product is for in vitro research use only