

CRK Monoclonal Antibody

Catalog No: #27127

Package Size: #27127-1 50ul #27127-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


| | |
|-----------------------|---|
| Product Name | CRK Monoclonal Antibody |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 3H7-E5-H8 |
| Isotype | IgG2b |
| Purification | Affinity purified |
| Applications | ELISA WB |
| Species Reactivity | Hu D Rb Ms |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Purified recombinant human CRK protein fragments expressed in E.coli. |
| Target Name | CRK |
| Other Names | p38, CRKII |
| Accession No. | Swiss-Prot:P05362Gene ID:3383 |
| Uniprot | P05362 |
| GeneID | 3383; |
| Formulation | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol |

Application Details

Predicted MW: 34kd

Western blotting: 1:1000

Images



Western blot detection of CRK in LNCAP,A431

Background

The Crk-I and Crk-II forms differ in their biological activities. Crk-II has less transforming activity than Crk-I. Crk-II mediates attachment-induced

MAPK8 activation, membrane ruffling and cell motility in a Rac-dependent manner.

Note: This product is for in vitro research use only