

## UCHL1 Monoclonal Antibody

Catalog No: #27128

Package Size: #27128-1 50ul #27128-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	UCHL1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7B8-G5-B8
Isotype	IgG2b
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human UCHL1 protein fragments expressed in E.coli.
Target Name	UCHL1
Other Names	PARK5, PGP95, PGP9.5, Uch-L1, PGP 9.5
Accession No.	Swiss-Prot:P24864 Gene ID:898
Uniprot	P24864
GeneID	898;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

## Application Details

Predicted MW: 25kd

Western blotting: 1:1000

## Images



Western blot detection of UCHL1 in U251,IMR32,Raji&COS7 cell lysates and using UCHL1 antibody (1:1000 diluted).  
 Predicted band size: 25KDa  
 Observed band size: 28KDa

## Background

Ubiquitin-protein hydrolase involved both in the processing of ubiquitin precursors and of ubiquitinated proteins. This enzyme is a thiol protease that

recognizes and hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin.

---

Note: This product is for in vitro research use only