

IL-5 Antibody FITC Conjugated

Catalog No: #C02398F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-5 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 70-120 132 derived from mouse IL-5
Conjugates	FITC
Target Name	IL-5
Other Names	IL-5; Interleukin-5; B-cell growth factor II; BCGF-II; Cytotoxic T-lymphocyte inducer; Eosinophil differentiation factor; T-cell replacing factor; TRF; IL5
Accession No.	Swiss-Prot#P04401NCBI Gene ID16191
Uniprot	P04401
GeneID	16191;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Factor that induces terminal differentiation of late-developing B-cells to immunoglobulin secreting cells.

Note: This product is for in vitro research use only