

# Recombinant Human Granulocyte Colony Stimulating Factor (rHu G-CSF)

Catalog No: #70202

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Granulocyte Colony Stimulating Factor (rHu G-CSF)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Species Reactivity	Hu
Target Name	rHu G-CSF
Other Names	IL2, T-cell Growth Factor, TCGF, Aldesleukin
Accession No.	accession:P60568 GeneID:3558
Uniprot	P60568
GeneID	3558;
Calculated MW	Approximately 15.4 kDa, a sing
SDS-PAGE MW	Sterile Filtered White lyophil
Target Sequence	APTSSSTKKT QLQLEHLLLD LQMILNGINN YKNPKLTRML TFKFYMPKKA TELKHLQCLE EELKPLEEVL NLAQSKNFHL RPRDLISNIN VIVLELKGSE TTFMCEYADE TATIVEFLNR WITFCQSIIS TLT
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 5.0.
Storage	but should be kept at -20 °C for long term storage

## Background

preferably desiccated. Upon reconstitution

## References

the preparation is stable for up to one week at 2-8 °C. For maximal stability

Note: This product is for in vitro research use only