

Recombinant Human BMP-2

Catalog No: #70802

Package Size: #70802-1 10ug #70802-2 100ug #70802-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BMP-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Species Reactivity	Hu
Target Name	rHu BMP-2
Other Names	BMP-2A
Accession No.	accession:P12643 GeneID:650
Uniprot	P12643
GeneID	650;
Calculated MW	Approximately 26.0 kDa, a homo
SDS-PAGE MW	Sterile Filtered White lyophil
Target Sequence	MQAKHKQRKR LKSSCKRHPL YVDFSDVGWN DWIVAPPGYH AFYCHGECPF PLADHLNSTN HAIVQTLVNS VNSKIPKACC VPTELSAISM LYLDENEKVV LKNYQDMVVE GCGCR
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution containing 10 mM sodium citrate pH 3.5.
Storage	but should be kept at -20 °C for long term storage

Background

preferably desiccated. Upon reconstitution

References

the preparation is stable for up to one week at 2-8 °C. For maximal stability

Note: This product is for in vitro research use only