

Rabbit anti-Human IgG secondary antibody HRP-conjugated

Catalog No: #L3052



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabbit anti-Human IgG secondary antibody HRP-conjugated
Brief Description	Secondary Antibody
Clonality	Polyclonal
Applications	ELISA WB IHC
Species Reactivity	Hu
Immunogen Description	The antiserum was developed in rabbit using the human IgG as the immunogen.
Conjugates	HRP
Target Name	Human IgG
Formulation	HRP-labeled Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 50% glycerol.
Storage	Store at -20 °C/1 year

Application Details

ELISA: 1:5,000 - 1:20,000 WB: 1:5,000 - 1:20,000 IHC: 1:500 - 1:1,000

Note: This product is for in vitro research use only