

ZBED2 Antibody FITC Conjugated

Catalog No: #C02668F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZBED2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide derived from human ZBED2
Conjugates	FITC
Target Name	ZBED2
Other Names	MGC10796; ZBED2; ZBED2_HUMAN; Zinc finger BED domain containing protein 2; Zinc finger BED domain-containing protein 2; Zinc finger, BED domain containing 2.
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

ZBED2

Note: This product is for in vitro research use only