

MAP2 MAP-2a.b.c Antibody FITC Conjugated

Catalog No: #C02791F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAP2 MAP-2a.b.c Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human MAP2
Conjugates	FITC
Target Name	MAP2 MAP-2a.b.c
Other Names	MAP2A; MAP2B; MAP2C; Microtubule-associated protein 2; MAP-2; MAP2
Accession No.	Swiss-Prot#P11137NCBI Gene ID4133
Uniprot	P11137
GeneID	4133;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

The exact function of MAP2 is unknown but MAPs may stabilize the microtubules against depolymerization. They also seem to have a stiffening effect on microtubules.

Note: This product is for in vitro research use only