

ADAMTS12 Antibody HRP Conjugated

Catalog No: #C02851H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADAMTS12 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 310-360 1594 derived from human ADAMTS12
Conjugates	HRP
Target Name	ADAMTS12
Other Names	PRO4389; A disintegrin and metalloproteinase with thrombospondin motifs 12; ADAM-TS 12; ADAM-TS12; ADAMTS-12; ADAMTS12; UNQ1918 PRO4389
Accession No.	Swiss-Prot#P58397NCBI Gene ID81792
Uniprot	P58397
GeneID	81792;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Metalloprotease that may play a role in the degradation of COMP. Cleaves also alpha-2 macroglobulin and aggrecan. Has anti-tumorigenic properties.

Note: This product is for in vitro research use only