

IRAK2 Antibody HRP Conjugated

Catalog No: #C02897H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IRAK2 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human IRAK2
Conjugates	HRP
Target Name	IRAK2
Other Names	Interleukin-1 receptor-associated kinase-like 2; IRAK-2; IRAK2
Accession No.	NCBI Gene ID3656
Uniprot	O43187
GeneID	3656;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

IRAK or Interleukin-1 Receptor-associated Kinase 1, is one of two putative serine threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. This protein is partially responsible for IL-1-induced upregulation of the ubiquitous transcription factor NF-kappa B.

Note: This product is for in vitro research use only