

## RASD1 Antibody HRP Conjugated

Catalog No: #C02898H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RASD1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 30-60 281 derived from human RASD1
Conjugates	HRP
Target Name	RASD1
Other Names	AGS1; DEXRAS1; MGC:26290; Dexamethasone-induced Ras-related protein 1; Activator of G-protein signaling 1; RASD1
Accession No.	Swiss-Prot#Q9Y272NCBI Gene ID51655
Uniprot	Q9Y272
GeneID	51655;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

## Background

Small GTPase. Negatively regulates the transcription regulation activity of the APBB1 FE65-APP complex via its interaction with APBB1 FE65 (By similarity).

Note: This product is for in vitro research use only