

## Gastroke 1 Antibody FITC Conjugated

Catalog No: #C02920F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gastroke 1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 130-160 199 derived from human Gastroke 1
Conjugates	FITC
Target Name	Gastroke 1
Other Names	FOV; CA11; AMP18; BRICD1; foveolin; Gastroke-1; 18 kDa antrum mucosa protein; AMP-18; Protein CA11; GKN1; UNQ489 PRO1005
Accession No.	Swiss-Prot#Q9NS71NCBI Gene ID56287
Uniprot	Q9NS71
GeneID	56287;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

## Background

Has mitogenic activity and may be involved in maintaining the integrity of the gastric mucosal epithelium.

Note: This product is for in vitro research use only