

C7 Antibody HRP Conjugated

Catalog No: #C03221H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---------------------------------------------------------------------------------|
| Product Name | C7 Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | WB IHC-P IHC-F ICC |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide derived from human C7 |
| Conjugates | HRP |
| Target Name | C7 |
| Other Names | C7; CO7_HUMAN; complement component 7; Complement component C7. |
| Excitation Emission | N A |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

C7

Note: This product is for in vitro research use only