

IDO Antibody FITC Conjugated

Catalog No: #C03438F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IDO Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 110-160 403 derived from human IDO
Conjugates	FITC
Target Name	IDO
Other Names	IDO; INDO; IDO-1; Indoleamine 2,3-dioxygenase 1; Indoleamine-pyrrole 2,3-dioxygenase; IDO1
Accession No.	Swiss-Prot#P14902NCBI Gene ID3620
Uniprot	P14902
GeneID	3620;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

Catalyzes the cleavage of the pyrrol ring of tryptophan and incorporates both atoms of a molecule of oxygen.

Note: This product is for in vitro research use only