

AQP5 Antibody FITC Conjugated

Catalog No: #C03486F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AQP5 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 227-265 265 derived from mouse AQP5
Conjugates	FITC
Target Name	Aquaporin 5
Other Names	Aquaporin-5; AQP-5; Aqp5
Accession No.	Swiss-Prot#Q9WTY4NCBI Gene ID11830
Uniprot	Q9WTY4
GeneID	11830;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Forms a water-specific channel. Implicated in the generation of saliva, tears, and pulmonary secretions.

Note: This product is for in vitro research use only