

K19 CYFRA21-1 Antibody FITC Conjugated

Catalog No: #C03521F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	K19 CYFRA21-1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 310-360 400 derived from human K19 CYFRA21-1
Conjugates	FITC
Target Name	K19 CYFRA21-1
Other Names	K19; CK19; K1CS; Keratin, type I cytoskeletal 19; Cytokeratin-19; CK-19; Keratin-19; KRT19
Accession No.	Swiss-Prot#P08727NCBI Gene ID3880
Uniprot	P08727
GeneID	3880;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Involved in the organization of myofibers. Together with KRT8, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.

Note: This product is for in vitro research use only