

K19 CYFRA21-1 Antibody HRP Conjugated

Catalog No: #C03521H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | K19 CYFRA21-1 Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | WB IHC-P IHC-F |
| Species Reactivity | Hu |
| Immunogen Description | KLH conjugated synthetic peptide aa 310-360 400 derived from human K19 CYFRA21-1 |
| Conjugates | HRP |
| Target Name | K19 CYFRA21-1 |
| Other Names | K19; CK19; K1CS; Keratin, type I cytoskeletal 19; Cytokeratin-19; CK-19; Keratin-19; KRT19 |
| Accession No. | Swiss-Prot#P08727NCBI Gene ID3880 |
| Uniprot | P08727 |
| GeneID | 3880; |
| Excitation Emission | N A |
| Cell Localization | Secreted |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Involved in the organization of myofibers. Together with KRT8, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.

Note: This product is for in vitro research use only