

GLP-1R Antibody HRP Conjugated

Catalog No: #C03524H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | GLP-1R Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | WB IHC-P IHC-F ICC |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 115-145 463 derived from rat GLP-1R |
| Conjugates | HRP |
| Target Name | GLP-1R |
| Other Names | GLP-1R; GLP1Rc; Glucagon-like peptide 1 receptor; GLP-1 receptor; GLP-1-R; Glp1r |
| Accession No. | Swiss-Prot#O35659NCBI Gene ID14652 |
| Uniprot | O35659 |
| GeneID | 14652; |
| Excitation Emission | N A |
| Cell Localization | Extracellular, Helical |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

This is a receptor for glucagon-like peptide 1. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Note: This product is for in vitro research use only