

## ENPP3 CD203c Antibody FITC Conjugated

Catalog No: #C03531F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ENPP3 CD203c Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human ENPP3: 40-140/875
Conjugates	FITC
Target Name	ENPP3 CD203c
Other Names	CD23c; A1876438; Ectonucleotide pyrophosphatase phosphodiesterase family member 3; E-NPP 3; Phosphodiesterase I beta; PD-Ibeta; Phosphodiesterase I nucleotide pyrophosphatase 3; Enpp3
Accession No.	Swiss-Prot#Q6DYE8NCBI Gene ID209558
Uniprot	Q6DYE8
GeneID	209558;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt 2ug /test IF 1:200-800

## Background

Cleaves a variety of phosphodiester and phosphosulfate bonds including deoxynucleotides, nucleotide sugars, and NAD.

Note: This product is for in vitro research use only