

## CK20 Antibody FITC Conjugated

Catalog No: #C03541F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CK20 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 180-225 424 derived from human CK20
Conjugates	FITC
Target Name	CK20
Other Names	K20; CD20; CK20; CK-20; KRT21; Keratin, type I cytoskeletal 20; Cytokeratin-20; Keratin-20; Protein IT; KRT20
Accession No.	Swiss-Prot#P35900NCBI Gene ID54474
Uniprot	P35900
GeneID	54474;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Plays a significant role in maintaining keratin filament organization in intestinal epithelia. When phosphorylated, plays a role in the secretion of mucin in the small intestine (By similarity).

Note: This product is for in vitro research use only