

FBXO48 Antibody HRP Conjugated

Catalog No: #C03552H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FBXO48 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 80-130 155 derived from human FBXO48
Conjugates	HRP
Target Name	FBXO48
Other Names	F-box only protein 48; FBXO48; FBX48
Accession No.	Swiss-Prot#Q5FWF7NCBI Gene ID554251
Uniprot	Q5FWF7
GeneID	554251;
Excitation Emission	N A
Cell Localization	Intracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

F-box containing protein.

Note: This product is for in vitro research use only