

CD15 Fut4 SSEA-1 Antibody FITC Conjugated

Catalog No: #C03638F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD15 Fut4 SSEA-1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 250-300 433 derived from mouse CD15
Conjugates	FITC
Target Name	CD15 Fut4 SSEA-1
Other Names	FAL; LeX; CD15; Ssea1; SSEA-1; FucT-IV; AI451562; Alpha-(1,3)-fucosyltransferase 4; Fucosyltransferase 4; Fucosyltransferase IV; Fuc-TIV; Galactoside 3-L-fucosyltransferase; Fut4; Elft
Accession No.	Swiss-Prot#Q11127NCBI Gene ID14345
Uniprot	Q11127
GeneID	14345;
Excitation Emission	494nm 518nm
Cell Localization	Intracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

May catalyze alpha-1,3 glycosidic linkages involved in the expression of Lewis X SSEA-1 and VIM-2 antigens.

Note: This product is for in vitro research use only