

PAI-1 Antibody FITC Conjugated

Catalog No: #C03640F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PAI-1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 330-380 402 derived from human PAI-1
Conjugates	FITC
Target Name	PAI-1
Other Names	PAI; PAI1; PAI-1; PLANH1; Plasminogen activator inhibitor 1; Endothelial plasminogen activator inhibitor; Serpin E1; SERPINE1
Accession No.	Swiss-Prot#P05121NCBI Gene ID5054
Uniprot	P05121
GeneID	5054;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Serine protease inhibitor. This inhibitor acts as 'bait' for tissue plasminogen activator, urokinase, protein C and matriptase-3 TMPRSS7. Its rapid interaction with PLAT may function as a major control point in the regulation of fibrinolysis.

Note: This product is for in vitro research use only