

## TRIM62 Antibody HRP Conjugated

Catalog No: #C03643H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TRIM62 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 150-200 475 derived from human TRIM62
Conjugates	HRP
Target Name	TRIM62
Other Names	DEAR1; E3 ubiquitin-protein ligase TRIM62; Tripartite motif-containing protein 62; TRIM62
Accession No.	Swiss-Prot#Q9BVG3NCBI Gene ID55223
Uniprot	Q9BVG3
GeneID	55223;
Excitation Emission	N A
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

## Background

E3 ubiquitin ligase whose activity is dependent on E2 ubiquitin-conjugating enzyme UBE2D2.

Note: This product is for in vitro research use only