

## HPV16 E6 + HPV18 E6 Antibody HRP Conjugated

Catalog No: #C03653H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HPV16 E6 + HPV18 E6 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu
Crossing Reactivity	HPV16 HPV18
Immunogen Description	KLH conjugated synthetic peptides from HPV16 E6 + HPV18 E6
Conjugates	HRP
Target Name	HPV16 E6 + HPV18 E6
Other Names	Protein E6; E6
Accession No.	Swiss-Prot#P03126, P06463NCBI Gene ID1489078, 1489088
Uniprot	P03126
GeneID	1489078;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200

## Background

Human Papilloma virus (HPV) is responsible for the development of type and tissue specific papilloma in the oral cavity and in the larynx. In the skin papilloma virus can cause different types of warts and in the ano-genital region it is associated with condyloma or carcinoma. HPV types 16 and 18 seem to be responsible for the development of some carcinoma like cervical carcinoma.

Note: This product is for in vitro research use only