

ABCD1 CCL22 Antibody FITC Conjugated

Catalog No: #C03682F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABCD1 CCL22 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 60-90 93 derived from mouse ABCD1
Conjugates	FITC
Target Name	ABCD1 CCL22
Other Names	MDC; DCBCK; ABCD-1; Scya22; C-C motif chemokine 22; Activated B and dendritic cell-derived; CC chemokine ABCD-1; Small-inducible cytokine A22; Ccl22; Abcd1
Accession No.	Swiss-Prot#O88430NCBI Gene ID20299
Uniprot	O88430
GeneID	20299;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Chemotactic for activated T-lymphocytes. May play an important role in the collaboration of dendritic cells and B-lymphocytes with T-cells in immune responses.

Note: This product is for in vitro research use only