

STAT6 (Phospho-Tyr641) Antibody HRP Conjugated

Catalog No: #C03703H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	STAT6 (Phospho-Tyr641) Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic phosphopeptide aa 630-645 847 derived from human STAT6 around the phosphorylation site of tyrosine Tyr641
Conjugates	HRP
Target Name	STAT6 Tyr641
Other Names	STAT6B; STAT6C; D12S1644; IL-4-STAT; Signal transducer and activator of transcription 6; IL-4 Stat; STAT6
Accession No.	Swiss-Prot#P42226NCBI Gene ID6778
Uniprot	P42226
GeneID	6778;
Excitation Emission	N A
Cell Localization	Cytoplasm, Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Carries out a dual function: signal transduction and activation of transcription. Involved in IL4 interleukin-4- and IL3 interleukin-3-mediated signaling.

Note: This product is for in vitro research use only