

IL-7 Antibody FITC Conjugated

Catalog No: #C03716F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-7 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 50-100 154 derived from mouse IL-7
Conjugates	FITC
Target Name	IL-7
Other Names	IL-7; hlb368; A630026I06Rik; Interleukin-7; IL7
Accession No.	Swiss-Prot#P10168NCBI Gene ID16196
Uniprot	P10168
GeneID	16196;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Hematopoietic growth factor capable of stimulating the proliferation of lymphoid progenitors. It is important for proliferation during certain stages of B-cell maturation.

Note: This product is for in vitro research use only