

Smad1 (Ab-206) Antibody

Catalog No: #21579

Package Size: #21579-1 50ul #21579-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Smad1 (Ab-206) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Target Name	Smad1
Other Names	BSP1, JV41, MADH1, MADR1,
Accession No.	Swiss-Prot:Q15797Gene ID:4086
Uniprot	Q15797
GeneID	4086;
Concentration	1.0mg/ml

Application Details

Predicted MW: 60kd

Western blotting: 1: 500~1: 1000

Images



Western blot analysis of extracts from HeLa and 3T3 cells using Smad1 (Ab-206) Antibody #21579.

Note: This product is for in vitro research use only