

## STEAP1 Antibody FITC Conjugated

Catalog No: #C03783F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	STEAP1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 290-339 339 derived from human STEAP1
Conjugates	FITC
Target Name	STEAP1
Other Names	STEAP; PRSS24; Metalloreductase STEAP1; Six-transmembrane epithelial antigen of prostate 1; STEAP1
Accession No.	Swiss-Prot#Q9UHE8NCBI Gene ID26872
Uniprot	Q9UHE8
GeneID	26872;
Excitation Emission	494nm 518nm
Cell Localization	Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Metalloreductase that has the ability to reduce both Fe(3+) to Fe(2+) and Cu(2+) to Cu(1+). Uses NAD(+) as acceptor (By similarity).

Note: This product is for in vitro research use only