

PLAUR Antibody FITC Conjugated

Catalog No: #C03800F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PLAUR Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 220-270 335 derived from human PLAUR
Conjugates	FITC
Target Name	PLAUR
Other Names	CD87; UPAR; URKR; U-PAR; Urokinase plasminogen activator surface receptor; Monocyte activation antigen Mo3; PLAUR; MO3
Accession No.	Swiss-Prot#Q03405NCBI Gene ID5329
Uniprot	Q03405
GeneID	5329;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Acts as a receptor for urokinase plasminogen activator. Plays a role in localizing and promoting plasmin formation. Mediates the proteolysis-independent signal transduction activation effects of U-PA. It is subject to negative-feedback regulation by U-PA which cleaves it into an inactive form.

Note: This product is for in vitro research use only