

RFP-Tag Mouse Monoclonal Antibody

Catalog No: #T608

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RFP-Tag Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Affinity-chromatography
Applications	WB
Specificity	RFP-Tag mouse monoclonal antibody detects endogenous levels of internal, C-terminal, and N-terminal RFP recombinant proteins.
Immunogen Type	Recombinant Protein
Immunogen Description	Full length RFP protein.
Target Name	RFP-Tag
SDS-PAGE MW	N/A
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C. Stable for 12 months from date of receipt

Application Details

WB: 1:3000~1:10000

Images



1ug RFP fusion protein+ Primary antibody dilution at
1γ • 1:3,000 2γ • 1:5,000 3γ • 1:10,000

Note: This product is for in vitro research use only