

REV3L Antibody HRP Conjugated

Catalog No: #C03852H

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	REV3L Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F ICC
Species Reactivity	Hu Ms
Immunogen Description	KLH conjugated synthetic peptide aa 2750-2790 3130 derived from human REV3L
Conjugates	HRP
Target Name	REV3L
Other Names	POLZ; REV3; DNA polymerase zeta catalytic subunit; Protein reversionless 3-like; REV3-like; hREV3; REV3L
Accession No.	Swiss-Prot#O60673NCBI Gene ID5980
Uniprot	O60673
GeneID	5980;
Excitation Emission	N A
Cell Localization	Cytoplasm, Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

Interacts with MAD2L2 to form the error prone DNA polymerase zeta involved in translesion DNA synthesis.

Note: This product is for in vitro research use only