

E Tag Antibody HRP Conjugated

Catalog No: #C03957H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E Tag Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB
Species Reactivity	Recombinant fusion protein E-tagged
Crossing Reactivity	Recombinant fusion protein, E tag
Immunogen Description	KLH conjugated E Tag synthetic peptide (GAPVPYDPLEPR)
Conjugates	HRP
Target Name	E Tag
Other Names	GAPVPYDPLEPR tag antibody; GAPVPYDPLEPR epitope tag antibody; E epitope tag antibody.
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:10000-20000

Background

E Tag is a commonly used epitope tag engineered onto the N- or C- terminus of a protein of interest so that the tagged protein can be analyzed and visualized. Because of the small size of the epitope, it is unlikely to alter the activity of the cloned sequence.

Note: This product is for in vitro research use only