

Angiotensin II type 1A receptor Antibody FITC Conjugated

Catalog No: #C03965F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin II type 1A receptor Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 163-192 359 derived from human of human Angiotensin II type 1A receptor
Conjugates	FITC
Target Name	Angiotensin II type 1A receptor
Other Names	AT1; AG2S; AT1B; AT1R; AT1AR; AT1BR; AT2R1; HAT1R; AGTR1B; Type-1 angiotensin II receptor; Angiotensin II type-1 receptor; AGTR1; AGTR1A; AT2R1B
Accession No.	Swiss-Prot#P30556NCBI Gene ID185
Uniprot	P30556
GeneID	185;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Background

Receptor for angiotensin II. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system.

Note: This product is for in vitro research use only