

## Glypican 4 Antibody FITC Conjugated

Catalog No: #C03974F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glypican 4 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Glypican 4
Conjugates	FITC
Target Name	Glypican 4
Other Names	dJ900E8.1 glypican 4; Glypican proteoglycan 4; Gpc4; K Glypican; Kglypican; OTTHUMP00000024052; Secreted glypican 4; GPC4_HUMAN.
Accession No.	NCBI Gene ID2239
Uniprot	O75487
GeneID	2239;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

## Background

Glypicans are glycosylphosphatidylinositol-anchored heparan sulphate proteoglycans, which are thought to modulate signalling by growth factors, including fibroblast growth factors and Wnts. Expression of glypican-4 has been reported in the mouse brain at embryonic day 7 and later stages.

Note: This product is for in vitro research use only