

## GAPDH Antibody HRP Conjugated

Catalog No: #C03990H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | GAPDH Antibody HRP Conjugated   |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Isotype               | IgG   |
| Purification          | Purified by Protein A.  |
| Applications          | WB IHC-P IHC-F  |
| Species Reactivity    | Hu Ms Rt Rab  |
| Immunogen Description | GAPDH protein of rabbit   |
| Conjugates            | HRP   |
| Target Name           | GAPDH   |
| Other Names           | G3PD; GAPD; HEL-S-162eP; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH; CDABP47; OK SW-cl.12 |
| Accession No.         | Swiss-Prot#P04406NCBI Gene ID2597   |
| Uniprot               | P04406  |
| GeneID                | 2597;   |
| Excitation Emission   | N A   |
| Cell Localization     | Cytoplasm   |
| Concentration         | 1mg ml  |
| Formulation           | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.  |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.   |

## Application Details

WB=1:10000-20000 IHC-P=1:50-200 IHC-F=1:50-200

## Background

Has both glyceraldehyde-3-phosphate dehydrogenase and nitrosylase activities, thereby playing a role in glycolysis and nuclear functions, respectively. Participates in nuclear events including transcription, RNA transport, DNA replication and apoptosis. Nuclear functions are probably due to the nitrosylase activity that mediates cysteine S-nitrosylation of nuclear target proteins such as SIRT1, HDAC2 and PRKDC. Modulates the organization and assembly of the cytoskeleton. Facilitates the CHP1-dependent microtubule and membrane associations through its ability to stimulate the binding of CHP1 to microtubules (By similarity). Glyceraldehyde-3-phosphate dehydrogenase is a key enzyme in glycolysis that catalyzes the first step of the pathway by converting D-glyceraldehyde 3-phosphate (G3P) into 3-phospho-D-glyceroyl phosphate. Component of the GAIT (gamma interferon-activated inhibitor of translation) complex which mediates interferon-gamma-induced transcript-selective translation inhibition in inflammation processes. Upon interferon-gamma treatment assembles into the GAIT complex which binds to stem loop-containing GAIT elements in the 3'-UTR of diverse inflammatory mRNAs (such as ceruplasmin) and suppresses their translation.

Note: This product is for in vitro research use only