

LDL Antibody HRP Conjugated

Catalog No: #C03996H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LDL Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms
Immunogen Description	Full length plasma protein (Human)
Conjugates	HRP
Target Name	LDL
Other Names	AD2; LPG; APO-E; LDLQ5; Apolipoprotein E; APOE
Accession No.	Swiss-Prot#P02649NCBI Gene ID348
Uniprot	P02649
GeneID	348;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Mediates the binding, internalization, and catabolism of lipoprotein particles. It can serve as a ligand for the LDL (apo B E) receptor and for the specific apo-E receptor (chylomicron remnant) of hepatic tissues.

Note: This product is for in vitro research use only