

## Pectinesterase inhibitor 18 Antibody HRP Conjugated

Catalog No: #C04006H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pectinesterase inhibitor 18 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Arabidopsis Thaliana
Crossing Reactivity	Arabidopsis thaliana,
Immunogen Description	KLH conjugated synthetic peptide derived from Arabidopsis thaliana pectinesterase inhibitor 18 Pectinesterase
Conjugates	HRP
Target Name	Pectinesterase inhibitor 18
Other Names	Pectinesterase pectinesterase inhibitor 18; AtPMEpcrA; Contains: Pectinesterase inhibitor 18; Pectin methylesterase inhibitor 18; Contains:Bunctional pectinesterase 18 rRNA N-glycosylase; PE 18; Pectin methylesterase 18; Pectin methylesterase 4; AtPME4; PME18_ARATH.
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

## Background

Pectinesterase inhibitor 18

Note: This product is for in vitro research use only