

CNR2 CB2 Antibody FITC Conjugated

Catalog No: #C04063F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CNR2 CB2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 298-330 360 derived from human CNR2
Conjugates	FITC
Target Name	CNR2 CB2
Other Names	CB2; CX5; CB-2; Cannabinoid receptor 2; hCB2; CNR2; CB2A; CB2B
Accession No.	Swiss-Prot#P34972NCBI Gene ID1269
Uniprot	P34972
GeneID	1269;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Heterotrimeric G protein-coupled receptor for endocannabinoid 2-arachidonoylglycerol mediating inhibition of adenylate cyclase. May function in inflammatory response, nociceptive transmission and bone homeostasis.

Note: This product is for in vitro research use only