

Endothelial lipase Antibody FITC Conjugated

Catalog No: #C04075F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Endothelial lipase Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | IF |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 310-360 500 derived from human Endothelial lipase |
| Conjugates | FITC |
| Target Name | Endothelial lipase |
| Other Names | EL; EDL; PRO719; Endothelial lipase; Endothelial cell-derived lipase; LIPG; UNQ387 PRO719 |
| Accession No. | Swiss-Prot#Q9Y5X9NCBI Gene ID9388 |
| Uniprot | Q9Y5X9 |
| GeneID | 9388; |
| Excitation Emission | 494nm 518nm |
| Cell Localization | Secreted |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Background

Has phospholipase and triglyceride lipase activities. Hydrolyzes high density lipoproteins (HDL) more efficiently than other lipoproteins. Binds heparin.

Note: This product is for in vitro research use only