

## BAFF Antibody HRP Conjugated

Catalog No: #C04100H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	BAFF Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 270-309 309 derived from mouse BAFF
Conjugates	HRP
Target Name	BAFF
Other Names	ApoL related ligand TALL 1; B cell Activating Factor; B lymphocyte stimulator ; BLYS; CD 257; CD257; CD257 antigen; delta BAFF; Dendritic cell derived TNF like molecule; DTL; DTL precursor; TALL 1; TALL1; THANK; TNF and APOL related leukocyte expressed ligand 1; TNF homolog that activates apoptosis;
Accession No.	Swiss-Prot#Q9WU72NCBI Gene ID24099
Uniprot	Q9WU72
GeneID	24099;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200

## Background

The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptors TNFRSF13B TACI, TNFRSF17 BCMA, and TNFRSF13C BAFFR. This cytokine is expressed in B cell lineage cells, and acts as a potent B cell activator. It has been also shown to play an important role in the proliferation and differentiation of B cells. [provided by RefSeq].

Note: This product is for in vitro research use only