

IL-13Ra2 Antibody FITC Conjugated

Catalog No: #C04123F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-13Ra2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 80-130 380 derived from human IL-13Ra2
Conjugates	FITC
Target Name	IL-13Ra2
Other Names	CT19; IL-13R; IL13BP; CD213A2; Interleukin-13 receptor subunit alpha-2; IL-13 receptor subunit alpha-2; IL-13R subunit alpha-2; IL-13R-alpha-2; IL-13RA2; Interleukin-13-binding protein; IL13RA2; IL13R
Accession No.	Swiss-Prot#Q14627NCBI Gene ID3598
Uniprot	Q14627
GeneID	3598;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Binds as a monomer with high affinity to interleukin-13 (IL13), but not to interleukin-4 (IL4).

Note: This product is for in vitro research use only