

TNFSF14 Antibody HRP Conjugated

Catalog No: #C04124H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNFSF14 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 145-195 240 derived from human TNFSF14
Conjugates	HRP
Target Name	TNFSF14
Other Names	LTg; TR2; CD258; HVEM-L; LIGHT; Tumor necrosis factor ligand superfamily member 14; Herpes virus entry mediator ligand; HVEM-L; Herpesvirus entry mediator ligand; TNFSF14; UNQ391 PRO726
Accession No.	Swiss-Prot#O43557NCBI Gene ID8740
Uniprot	O43557
GeneID	8740;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P IHC-F

Background

Cytokine that binds to TNFRSF3 LTBR. Binding to the decoy receptor TNFRSF6B modulates its effects. Activates NFκB, stimulates the proliferation of T-cells, and inhibits growth of the adenocarcinoma HT-29. Acts as a receptor for Herpes simplex virus.

Note: This product is for in vitro research use only