

Neuritin Antibody FITC Conjugated

Catalog No: #C04125F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuritin Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 70-116 142 derived from human Nrn1
Conjugates	FITC
Target Name	Neuritin
Other Names	NRN; dJ380B8.2; Neuritin; NRN1
Accession No.	Swiss-Prot#Q9NPD7NCBI Gene ID51299
Uniprot	Q9NPD7
GeneID	51299;
Excitation Emission	494nm 518nm
Cell Localization	Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Promotes neurite outgrowth and especially branching of neuritic processes in primary hippocampal and cortical cells.

Note: This product is for in vitro research use only