

NFκB-p65(Phospho-Ser468) Antibody

Catalog No: #11013

Package Size: #11013-1 50ul #11013-2 100ul


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NFκB-p65(Phospho-Ser468) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB IF
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Target Name	NFκB-p65
Modification	Phospho
Other Names	p65, NFKB3,
Accession No.	Swiss-Prot:Q04206Gene ID:5970
Uniprot	Q04206
GeneID	5970;
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use

Images



Immunofluorescence staining of methanol-fixed HeLa cells using NF-κB p65 (phospho-Ser468) antibody (#11013, Red).



Western blot analysis of extracts using NF-κB p65 (phospho-Ser468) antibody (#11013).

Note: This product is for in vitro research use only